THE TELEFLEX
ACADEMY


MC-XXXXXX

Focus: Le site IO humérus proximal

Stéphane Petitot | 10-05-2016

Teleflex®

SOMMAIRE:

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1. Internal study report. Protocol 2013-06: Clinical Studies to Determine the Optimal Technique to Identify the Proximal Humerus Intraosseous Vascular Access Insertion Site. Vidacare Corporation, May 2013.

2. Hoskins SL, Kramer GC, Stephens CT, Zachariah BS. Efficacy of epinephrine delivery via the intraosseous humeral head route during CPR. Circulation. 2006;(114):1204.

3. Hoskins SL, Zachariah BS, Copper N, Kramer GC. Comparison of intraosseous proximal humerus and sternal routes for drug delivery during CPR. Circulation. 2007;(116).

4. Hoskins SL, do Nascimento P, Lima RM, Espana-Tenorio JM, Kramer GC. Pharmacokinetics of intraosseous and central venous drug delivery during cardiopulmonary resuscitation. Resuscitation. 2012 Jan;83(1):107-12

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Annals of Medicine and Surgery (2015) 5:100–102

Contents lists available at ScienceDirect

Annals of Medicine and Surgery

Journal homepage: www.annalsjournal.com

Effects of humeral intraosseous versus intravenous epinephrine on pharmacokinetics and return of spontaneous circulation in a porcine cardiac arrest model: A randomized control trial

Don Johnson^{A, B, *}, Jose Garcia-Blanco^B, James Burgert^{A, B}, Lawrence Fulton^C, Patrick Kadlak^A, Katherine Perry^A, Jeffrey Burke^A

^ANorthwestern University, United States Army Graduate Program in Anesthetic Nursing, 3000 Regis #10, Suite 112, 1007 Van Housen, IL 60224, United States

^BSt. Charles Healthcare, 207 Madison Ave, Suite 402, Schaumburg, IL 60196, United States

^CState Illinois University, North College of Pharmacy, Dept. of Health Organization Management, Box 279, Lakeview, IL 60049, United States

HIGHLIGHTS

- No difference in concentration maximum (C_{max}) and time to maximum concentration (T_{max}) in epinephrine between humeral intraosseous and intravenous routes of administration over time.
- Humeral intraosseous delivers higher concentration than intravenous at 30 s after administration of epinephrine.
- Humeral intraosseous facilitates rapid delivery of epinephrine during cardiac arrest.
- Use of humeral intraosseous had higher number of subjects survived.

ARTICLE INFO

Article history:
Received 30 July 2015
Received in revised form 17 August 2015
Accepted 19 August 2015

Keywords:
Intraosseous
Return of spontaneous circulation
Epinephrine
Pharmacokinetics
Resuscitation

ABSTRACT

Cardiopulmonary resuscitation (CPR), defibrillation, and epinephrine administration are pillars of advanced cardiac life support (ACLS). Intraosseous (IO) access is an alternative route for epinephrine administration when intravenous (IV) access is unavailable. Previous studies indicate the pharmacokinetics of epinephrine administration via IO and IV routes differ, but it is not known if the difference influences return of spontaneous circulation (ROSC). The purpose of this prospective, experimental study was to determine the effects of humeral IO (HIO) and IV epinephrine administration during cardiac arrest on pharmacokinetics, ROSC, and odds of survival. Swine (n = 21) were randomized into 3 groups: humeral IO (HIO), peripheral IV (PIV) and CPR/defibrillation control. Cardiac arrest was induced under general anesthesia. The chest remained in arrest for 2 min without intervention. Chest compressions were initiated and continued for 2 min. Epinephrine was administered and arterial blood samples collected for pharmacokinetic analysis over 4 min. Defibrillation and epinephrine administration proceeded according to ACLS guidelines continuing for 20 min or until ROSC. Seven HIO swine, 6 PIV swine, and no control swine had ROSC. There were no significant differences in ROSC, maximum concentration, except at 30 s, and time-to-concentration-maximum between the HIO and IV groups. Significant differences existed between the experimental groups and the control. The HIO delivers a higher concentration of epinephrine than the IV route at 30 s which may be a survival advantage. Clinicians may consider using the IO route to administer epinephrine during CA when there is no preexisting IV access or when IV access is unavailable.

Published by Elsevier Ltd on behalf of JS Publishing Group Limited. This is an open access article under the CC BY-NC-ND license (<http://creativecommons.org/licenses/by-nc-nd/4.0/>).

1. Introduction

Incidence of death attributable to cardiovascular disease has declined over the past 15 years but still accounts for 1 of every 3

© 2015 Publishing services by Elsevier B.V. This is an open access article under the CC BY-NC-ND license (<http://creativecommons.org/licenses/by-nc-nd/4.0/>).

http://dx.doi.org/10.1016/j.amssurg.2015.08.010

2015-08-20 | Published by JS Publishing Group Limited. This is an open access article under the CC BY-NC-ND license (<http://creativecommons.org/licenses/by-nc-nd/4.0/>).

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Table 2
Cross-tabulation of group members by status.

		Groups			Total
		CPR with defib	HIO	IV	
Survive?	Yes	0	7	4	11
	No	7	0	3	10
Total		7	7	7	21


Fig. 1. Boxplots of Cmax by treatment group.


Fig. 2. Boxplots of Tmax by treatment group.

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
Fig. 3. Epinephrine concentration over time.

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

- Structures musculo tendineuse

Tendon bicipital

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Proximal Humerus Positioning / Anatomy Dissected (Cadaveric)

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

Fig. 1 Descriptive analysis of humeral IO success data.


1. D. Wampler, D. Schwartz, J Shumaker, S. Bolleter, R. Beckett, C. Manifold. Paramedics successfully perform humeral EZ-IO intraosseous access in adult out-of-hospital cardiac arrest patients, American Journal of Emergency Medicine (2011)

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
Des questions?

[Placeholder for additional 3rd party trademark owner]

Teleflex, the Teleflex logo, [List here in alphabetical order any other Teleflex trademarks found displayed in the document] are trademarks or registered trademarks of Teleflex Incorporated or its affiliates, in the U.S. and/or other countries.

© 20XX Teleflex Incorporated. All rights reserved. MC-XXXXXX

Merci

[Placeholder for additional 3rd party trademark owner]

Teleflex, the Teleflex logo, [List here in alphabetical order any other Teleflex trademarks found displayed in the document] are trademarks or registered trademarks of Teleflex Incorporated or its affiliates, in the U.S. and/or other countries.

© 20XX Teleflex Incorporated. All rights reserved. MC-XXXXXX